

Oxidation Stability Tester

Introduction:

Aluminium casting four hole block bath. outer body made of M.S with powder coated paints. Temperature control by digital temperature controller. Temperature range Ambient to 200°C. 4 Nos. Air flow meter for supply of oxygen in Oxidation tube.



590 OS PC 25-

Standard Working	Ambient to 200°C.
Oxidation Tube	4 Nos. Air flow meter.
Temperature controller	Microprocessor based PID controller.
Display	7 Segment LED Dual Display.
Insulation	High quality Glass or Ceramic blanket.

All electrical appliances work on 230VAC, 50Hz, single phase. Other Voltages available on request. Due to continuous improvement in design, the product supplied may have modified features.

Model	Temperature range	Controller
#595101MD	Ambient to 200°C.	Digital Temperature Controller

- 'S': If the external body is made of SS 304 (Stainless Steel 304), the model's name will end with S (e.g., Model: #xxxxxS).
- 'M': If the external body is made of mild steel, the model's name will end with M (e.g., Model: #xxxxxM).